

TECHNICAL DATA

Contamac's polymethyl methacrylate (PMMA) material is available in clear and yellow.

Our PMMA product has been on the market for many years with an excellent record of safety and consistency.

Please Note: Regulatory requirements and standards vary from country to country, and are constantly evolving. As a global company we want to be sure we provide you with detailed technical information, specific to your market, where appropriate, rather than using the condensed and simplified technical information on the website. If you need to use technical data for quality paperwork, or for a regulatory submission, please contact your account manager to obtain this precise and detailed information to support your regulatory requirements, we will be happy to help.

Material Characteristics

PROPERTY

PMMA

Water Absorption by Gravimetry (%)

≤1.5%

Refractive Index at 20°C, 589nm

1.49

10% Cut-Off Wavelength (nm)

377

Blue Light Absorber

Available

Please note: Some values may have been rounded for presentation purposes. Please contact your account manager for further details.